



GHygrocell

Humidity Test Chambers



Tests corrosion
To saturated humidity environment.

Corrosion tests by exposure to a saturated humidity environment

Standard Features:

- With Constant Temperature and Climate
- In accordance with the following rules:
ASTM D1735, ASTM D2247, BS 3900/F9, DIN 50017, DIN 55958, DIN 55991,
ISO11503, ISO 6270-2:2005, ISO 4541, ISO3231,
NFT 30-077 and other similar



Soluzioni per il Laboratorio e l'Industria
Laboratory and Industrial Solutions

Camere Simulazione Ambientale per la Ricerca, controllo Qualità e Produzione
Environmental Test Chambers for Research, Quality control and Production



- Structure made of corrosion-resistant material
- Crystal door for internal viewing of test pieces
- Heating by armoured resistors
- Excellent uniformity and thermal stability
- Structure made of corrosion-resistant material
- Crystal door for internal viewing of test pieces
- Heating by armoured resistors
- Excellent uniformity and thermal stability
- Device for lodging standard test pieces/sheets
- General switch bright green
- Protection fuses omnipolar
- Power Cable
- Safety plug with EEC interlock

Material supplied in equipment:

- N°01 User Manual
- N°01 Guarantee Certificate
- N°01 EC declaration



Soluzioni per il Laboratorio e l'Industria
Laboratory and Industrial Solutions

Camere Simulazione Ambientale per la Ricerca, controllo Qualità e Produzione
Environmental Test Chambers for Research, Quality control and Production