

GAutostill Distiller

Automatic Water Distiller

Automatic water distiller for the production of distilled water, in compliance with the European Pharmaceutical grade with conductivity below 2.3 μS at + 20 ° C; depending on the state / conditions of the water.

Operation

All the operational activities necessary for the distillation process are activated and controlled automatically, luminous lights indicate the various phases.

- consent to the supply of water,
- cooling water regulation,
- switching on the heating,
- switching off the heating in case of lack of water,
- shutdown of the heating and cooling water supply when the storage tank is full,
- ignition when the distillate has been taken.

Energy recovery system and reduction of water consumption, through heat exchange.



Technical features

Technical features	GAUTOSTILL07T
Hourly Production Liters (approx.)	7
Distillate conductivity μS / cm	2,3
Power supply V 50-60 Hz	400 3+N
Power W.	4800
Absorption A	21
Water inlet diameter mm	13
Cooling water outlet diameter mm	20
Distilled water outlet diameter mm	13
Dimensions LxWxH mm	710x300x500
Weight in Kg	20