

Doppio Frigorifero Professionale da Laboratorio



Ideale per la corretta conservazione a bassa temperatura.
 Ottimo isolamento termico su tutte le pareti.
 Interno in Acciaio Inox - Esterno verniciato.
 Ripiani estraibili e regolabili in altezza.
 Porta con guarnizioni magnetiche e chiave.
 Porta vetrata per l'osservazione interna.
 Microinterruttore Apertura ante.
 Pannello comandi con Sistema di Controllo della Temperatura Digitale con Display 0,1°C.
 Gestione degli Allarmi.
 Sonde di temperatura al Platino PT1000.

Compressore ermetico con Gas Ecologico.

Alimentazione 230V 1f 50Hz.

CODICE	GREFRY Doppio Frigorifero da Laboratorio	Dim. Interne LxPxH mm	Dim.Esterne LxPxH mm	Peso netto kg	Dim.Imballo LxPxH mm	Peso lordo kg
GREFRY-700VFF	GREFRY Doppio Vano 2x350L +2+12°C, Ventilato, 3+3 Ripiani, Digitale, All.AcusticoVisivo, DataLogger, USB	VS 580x665x720 VI 580x665x720	700x800x2040	183	750x850x2180	198



Professional Laboratory Double Door Refrigerator



Suitable for low temperature storage.
 Excellent thermal insulation on every wall.
 Stainless Steel internal chamber.
 Galvanized Steel exterior, PVC coated.
 Extractable and height adjustable shelves.
 Door with magnetic gaskets and lock with key.
 Glass door for internal observation.
 Microswitch at door opening.
 Control panel with digital temperature control system and 0,1°C display.
 Alarm control system.
 Platinum PT1000 temperature probes.

Hermetic compressor with ecofriendly gas.

Power supply 230V 1f 50Hz

CODE	GREFRY Laboratory Double Door Refrigerator	Interior Dim. WxDxH mm	Exterior Dim. WxDxH mm	Net Weight kg	Packaging Dim. WxDxH mm	Gross Weight kg
GREFRY-700VFF	GREFRY Double Compartment 2x350L +2+12°C, Ventilated, 3+3 Shelves, Digital, Visual Acoustic Alarm, DataLogger, USB	UC 580x665x720 LC 580x665x720	700x800x2040	183	750x850x2180	198

