

## Frigocongelatore Professionale da Laboratorio



Ideale per la corretta conservazione a bassa temperatura. Ottimo isolamento termico su tutte le pareti. Interno in Acciaio Inox - Esterno verniciato. Ripiani estraibili e regolabili in altezza. Porta con guarnizioni magnetiche e chiave. Versione con Porta vetrata per l'osservazione interna. Microinterruttore Apertura ante. Pannello comandi con Sistema di Controllo della Temperatura Digitale con Display 0,1°C. Gestione degli Allarmi. Sonde di temperatura al Platino PT1000. Compressore ermetico con Gas Ecologico.



Alimentazione 230V 1f 50Hz.

CODICE	REFRY Frigocongelatore da Laboratorio	Dim. Interne LxPxH mm	Dim.Esterne LxPxH mm	Peso netto kg	Dim.Imballo LxPxH mm	Peso lordo kg
<b>REFRY-700XCF</b>	REFRY Doppio Vano Frigorifero 350L +2+12°C, Ventilato, 4 Ripiani, Congelatore 350L -30-10°C, Digitale, All.AcusticoVisivo, DataLogger, USB	TN 580x670x685 BT 580x670x685	700x800x1990	180	750x850x2180	193
<b>REFRY-700VCF</b>	REFRY Doppio Vano Frigorifero 350L +2+12°C, Porta Vetrata, Ventilato, 4 Ripiani, Congelatore 350L -30-10°C, Digitale, All.AcusticoVisivo, DataLogger, USB	TN 580x670x685 BT 580x670x685	700x800x1990	183	750x850x2180	198



## Professional Laboratory Refrigerator and Freezer



Suitable for low temperature storage.  
 Excellent thermal insulation on every wall.  
 Stainless Steel internal chamber.  
 Galvanized Steel exterior, PVC coated.  
 Extractable and height adjustable shelves.  
 Door with magnetic gaskets and lock with key.  
 Available with glass door for internal observation.  
 Microswitch at door opening.  
 Control panel with digital temperature control system and 0,1°C display.  
 Alarm control system.  
 Platinum PT1000 temperature probes.  
 Hermetic compressor with ecofriendly gas.



Power supply 230V 1f 50Hz

CODE	GREFRY Laboratory Refrigerator and Freezer	Interior Dim. WxDxH mm	Exterior Dim. WxDxH mm	Net Weight kg	Packaging Dim. WxDxH mm	Gross Weight kg
<b>GREFRY-700XCF</b>	GREFRY Double Compartment Refrigerator 350L +2+12°C, Glass Door, Ventilated, 4 Shelves, Freezer 350L -30-10°C, Visual Acoustic Alarm, DataLogger, USB	NT 580x670x685 LT 580x670x685	700x800x1990	180	750x850x2180	193
<b>GREFRY-700VCF</b>	GREFRY Double Compartment Refrigerator 350L +2+12°C, Ventilated, 4 Shelves, Freezer 350L -30-10°C, Visual Acoustic Alarm, DataLogger, USB	NT 580x670x685 LT 580x670x685	700x800x1990	183	750x850x2180	198

