

## Frigorifero Professionale da Laboratorio



Ideale per la corretta conservazione a bassa temperatura.

Ottimo isolamento termico su tutte le pareti.

Interno in Acciaio Inox - Esterno verniciato.

Ripiani estraibili e regolabili in altezza.

Porta con guarnizioni magnetiche e chiave.

Versione con Porta vetrata per l'osservazione interna.

Microinterruttore Apertura ante.

Pannello comandi con Sistema di Controllo della Temperatura Digitale con Display 0,1°C.

Gestione degli Allarmi.

Sonde di temperatura al Platino PT1000.

Compressore ermetico con Gas Ecologico.

Alimentazione 230V 1f 50Hz.



CODICE	GREFRY Frigorifero da Laboratorio	Dim. Interne LxPxH mm	Dim.Esterne LxPxH mm	Peso netto kg	Dim.Imballo LxPxH mm	Peso lordo kg
<b>GREFRY-200XXX</b>	GREFRY 200L Frigorifero da Laboratorio +2+12°C, Ventilato, 3 Ripiani, Digitale, All.AcusticoVisivo, DataLogger, USB	480x465x920	600x600x1560	85	680x680x1790	96
<b>GREFRY-200VXX</b>	GREFRY 200L Frigorifero da Laboratorio +2+12°C, Porta Vetrata, Ventilato, 3 Ripiani, Digitale, All.AcusticoVisivo, DataLogger, USB	480x465x920	600x600x1560	91	680x680x1790	102

## Professional Laboratory Refrigerator



Suitable for low temperature storage.  
 Excellent thermal insulation on every wall.  
 Stainless Steel internal chamber.  
 Galvanized Steel exterior, PVC coated.  
 Extractable and height adjustable shelves.  
 Door with magnetic gaskets and lock with key.  
 Available with glass door for internal observation.  
 Microswitch at door opening.  
 Control panel with digital temperature control system and 0,1°C display.  
 Alarm control system.  
 Platinum PT1000 temperature probes.

Hermetic compressor with ecofriendly gas.

Power supply 230V 1f 50Hz



CODE	GREFRY Laboratory Refrigerator	Interior Dim. WxDxH mm	Exterior Dim. WxDxH mm	Net Weight kg	Packaging Dim. WxDxH mm	Gross Weight kg
<b>GREFRY-200XXX</b>	GREFRY 200L Laboratory Refrigerator +2+12°C, Ventilated, 3 Shelves, Digital, Visual Acoustic Alarm, DataLogger, USB	480x465x920	600x600x1560	85	680x680x1790	96
<b>GREFRY-200VXX</b>	GREFRY 200L Laboratory Refrigerator +2+12°C, Glass Door, Ventilated, 3 Shelves, Digital, Visual Acoustic Alarm, DataLogger, USB	480x465x920	600x600x1560	91	680x680x1790	102