

## Frigorifero Professionale da Laboratorio



Ideale per la corretta conservazione a bassa temperatura.  
 Ottimo isolamento termico su tutte le pareti.  
 Interno in Acciaio Inox - Esterno verniciato.  
 Ripiani estraibili e regolabili in altezza.  
 Porta con guarnizioni magnetiche e chiave.  
 Versione con Porta vetrata per l'osservazione interna.  
 Microinterruttore Apertura ante.  
 Pannello comandi con Sistema di Controllo della Temperatura Digitale con Display 0,1°C.  
 Gestione degli Allarmi.  
 Sonde di temperatura al Platino PT1000.

Compressore ermetico con Gas Ecologico.



Alimentazione 230V 1f 50Hz.

CODICE	GREFRY Frigorifero da Laboratorio	Dim. Interne LxPxH mm	Dim.Esterne LxPxH mm	Peso netto kg	Dim.Imballo LxPxH mm	Peso lordo kg
<b>GREFRY-300XXX</b>	GREFRY 300L Frigorifero da Laboratorio +2+12°C, Ventilato, 5 Ripiani, Digitale, All.AcusticoVisivo, DataLogger, USB	480x465x1320	600x600x1960	104	680x680x2150	117
<b>GREFRY-300VXX</b>	GREFRY 300L Frigorifero da Laboratorio +2+12°C, Porta Vetrata, Ventilato, 5 Ripiani, Digitale, All.AcusticoVisivo, DataLogger, USB	480x465x1320	600x600x1960	113	680x680x2150	126

## Professional Laboratory Refrigerator



Suitable for low temperature storage.  
 Excellent thermal insulation on every wall.  
 Stainless Steel internal chamber.  
 Galvanized Steel exterior, PVC coated.  
 Extractable and height adjustable shelves.  
 Door with magnetic gaskets and lock with key.  
 Available with glass door for internal observation.  
 Microswitch at door opening.  
 Control panel with digital temperature control system and 0,1°C display.  
 Alarm control system.  
 Platinum PT1000 temperature probes.  
 Hermetic compressor with ecofriendly gas.



Power supply 230V 1f 50Hz

CODE	GREFRY Laboratory Refrigerator	Interior Dim. WxDxH mm	Exterior Dim. WxDxH mm	Net Weight kg	Packaging Dim. WxDxH mm	Gross Weight kg
<b>GREFRY-300XXX</b>	GREFRY 300L Laboratory Refrigerator +2+12°C, Ventilated, 5 Shelves, Digital, Visual Acoustic Alarm, DataLogger, USB	480x465x1320	600x600x1960	104	680x680x2150	117
<b>GREFRY-300VXX</b>	GREFRY 300L Laboratory Refrigerator +2+12°C, Glass Door, Ventilated, 5 Shelves, Digital, Visual Acoustic Alarm, DataLogger, USB	480x465x1320	600x600x1960	113	680x680x2150	126

