

- Digital Control with Double Display
- Independent Safety Thermostat
- Temperature Range: from +10°C above ambient temperature to +260°C
- N°02 Height-Adjustable and Extractable Grid Shelves
  
- Internal dimensions 460x390x430mm WxDxH
- Volume 75 L
  
- External dimensions 600x570x730mm WxDxH
- Weight: 52 Kg
- Power: 1300 W Power Supply: 230V 1f 50Hz
  
- **State: Reconditioned with 6-month Warranty**
  
- Reference 21 s/n. 21.156

